

Product datasheet for **HP205475**

A RAF (ARAF) Human qPCR Primer Pair (NM_001654)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	369
Forward Sequence:	CCCACATTCCAAGTCACCAGCA
Reverse Sequence:	CCTCCCAGTAATAGCCTGAGTC
ACCN:	BC009028 , NM_001654 , NM_001654.1 , NM_001654.2 , NM_001654.3 , NM_001654.4 , BC002466 , BC002466.2 , BC007514 , BT019864 , NM_001654.5
UniProt ID:	P10398
Synonyms:	A-RAF; ARAF1; PKS2; RAFA1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)