

Product datasheet for **HP205481**

ARF4L (ARL4D) Human qPCR Primer Pair (NM_001661)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	379
Forward Sequence:	TGGACTCTGCTGGAAAGACCTC
Reverse Sequence:	CCACACTTGAAGGTGATGCCA
ACCN:	<u>NM_001661</u> , <u>NM_001661.1</u> , <u>NM_001661.2</u> , <u>NM_001661.3</u> , <u>BC000043</u> , <u>BC000043.2</u> , <u>BI547763</u> , <u>BM682844</u> , <u>BU729196</u> , <u>NM_001661.4</u>
UniProt ID:	<u>P49703</u>
Synonyms:	ARF4L; ARL6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)