

Product datasheet for **HP205483**

ARF6 Human qPCR Primer Pair (NM_001663)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	382
Forward Sequence:	CCAAGGTCTCATCTTCGTAGTGG
Reverse Sequence:	AGGCCTGCTTGTGGCGAAGA
ACCN:	<u>NM_001663</u> , <u>NM_001663.1</u> , <u>NM_001663.3</u> , <u>BC008918</u> , <u>BC008918.2</u> , <u>BC002952</u> , <u>BC030291</u> , <u>BM838197</u> , <u>BT019868</u> , <u>NM_001663.4</u>
UniProt ID:	<u>P62330</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)