

Product datasheet for **HP205491**

Asialoglycoprotein Receptor 1 (ASGR1) Human qPCR Primer Pair (NM_001671)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	432
Forward Sequence:	GAAGCAGTTCGTGTCTGACCTG
Reverse Sequence:	AGCGAGAGAACCAGTAGCAGCT
ACCN:	NM_001671 , NM_001671.1 , NM_001671.2 , NM_001671.3 , NM_001671.4 , BC032130 , BC032130.2 , BM818592 , BQ653541 , BT019873 , NM_001671.5
UniProt ID:	P07306
Synonyms:	ASGPR; ASGPR1; CLEC4H1; HL-1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)