

Product datasheet for **HP205504**

ATP6V0C Human qPCR Primer Pair (NM_001694)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	527
Forward Sequence:	GAATGACGACATCAGCCTCTACA
Reverse Sequence:	GATCATGCCACGAATAGTCGG
ACCN:	NM_001694 , NM_001694.1 , NM_001694.2 , NM_001694.3 , BC009290 , BC009290.2 , BC004537 , BC007389 , BC007759 , BI548787 , BT007155 , NM_001694.4
UniProt ID:	P27449
Synonyms:	ATP6C; ATP6L; ATPL; VATL; Vma3; VPPC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)