

Product datasheet for **HP205513**

bcl 6 (BCL6) Human qPCR Primer Pair (NM_001706)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	604
Forward Sequence:	CATGCAGAGATGTGCCTCCACA
Reverse Sequence:	TCAGAGAAGCGGCAGTCACACT
ACCN:	<u>NM_001706</u> , <u>NM_001706.1</u> , <u>NM_001706.2</u> , <u>NM_001706.3</u> , <u>NM_001706.4</u> , <u>BC142705</u> , <u>BC146796</u> , <u>BC150184</u> , <u>BX649185</u> , <u>NM_001706.5</u>
UniProt ID:	<u>P41182</u>
Synonyms:	BCL5; BCL6A; LAZ3; ZBTB27; ZNF51
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)