

Product datasheet for **HP205518**

Biglycan (BGN) Human qPCR Primer Pair (NM_001711)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	633
Forward Sequence:	TTGAACCTGGAGCCTTCGATGG
Reverse Sequence:	TTGGAGTAGCGAAGCAGGTCTT
ACCN:	<u>NM_001711</u> , <u>NM_001711.1</u> , <u>NM_001711.2</u> , <u>NM_001711.3</u> , <u>NM_001711.4</u> , <u>NM_001711.5</u> , <u>BC002416</u> , <u>BC004244</u> , <u>BQ180984</u> , <u>BT007323</u>
UniProt ID:	<u>P21810</u>
Synonyms:	DSPG1; MRLS; PG-S1; PGI; SEMDX; SLRR1A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)