

## Product datasheet for **HP205541**

### **C5R1 (C5AR1) Human qPCR Primer Pair (NM\_001736)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Gene ID:</b>	728
<b>Forward Sequence:</b>	ATCTGGTGCCAGAACTTCCGAG
<b>Reverse Sequence:</b>	CTGTAGTCCACGCCACACAACA
<b>ACCN:</b>	<u><a href="#">NM_001736</a></u> , <u><a href="#">NM_001736.1</a></u> , <u><a href="#">NM_001736.2</a></u> , <u><a href="#">NM_001736.3</a></u> , <u><a href="#">BC008982</a></u> , <u><a href="#">BC008982.1</a></u> , <u><a href="#">BE042476</a></u> , <u><a href="#">BI916506</a></u> , <u><a href="#">BM830753</a></u> , <u><a href="#">BT007358</a></u>
<b>UniProt ID:</b>	<u><a href="#">P21730</a></u>
<b>Synonyms:</b>	C5A; C5AR; C5R1; CD88
<b>Components:</b>	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix is stable for one year from date of shipping.



[View online »](#)