

Product datasheet for **HP205550**

Actin Regulatory Protein CAPG (CAPG) Human qPCR Primer Pair (NM_001747)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	822
Forward Sequence:	CAGGTGGAGATTGCTCACTGATGG
Reverse Sequence:	CTGGGCATTTGCCTTGTCAGCT
ACCN:	<u>NM_001747</u> , <u>NM_001747.1</u> , <u>NM_001747.2</u> , <u>NM_001747.3</u> , <u>BC000728</u> , <u>BC000728.2</u> , <u>BC014549</u> , <u>BG762114</u> , <u>BG762391</u> , <u>BU172582</u> , <u>BU626913</u>
UniProt ID:	<u>P40121</u>
Synonyms:	AFCP; HEL-S-66; MCP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)