

## Product datasheet for **HP205553**

### Caveolin 1 (CAV1) Human qPCR Primer Pair (NM\_001753)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	857
Forward Sequence:	CCAAGGAGATCGACCTGGTCAA
Reverse Sequence:	GCCGTCAAACCTGTGTGCCCT
ACCN:	<a href="#">NM_001753</a> , <a href="#">NM_001753.1</a> , <a href="#">NM_001753.2</a> , <a href="#">NM_001753.3</a> , <a href="#">NM_001753.4</a> , <a href="#">BC082246</a> , <a href="#">BC082246.1</a> , <a href="#">BC006432</a> , <a href="#">BC009685</a> , <a href="#">BG170208</a> , <a href="#">BT007143</a> , <a href="#">NM_001753.5</a>
UniProt ID:	<a href="#">Q03135</a>
Synonyms:	BSCL3; CGL3; LCCNS; MSTP085; PPH3; VIP21
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)