

Product datasheet for **HP205556**

CBR1 Human qPCR Primer Pair (NM_001757)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	873
Forward Sequence:	AGCTGGACATCGACGATCTGCA
Reverse Sequence:	TATGAAAGGGTGTGGGATCAGCA
ACCN:	NM_001757 , NM_001757.1 , NM_001757.2 , NM_001757.3 , BC015640 , BC015640.2 , BC002511 , BT019843 , NM_001757.4
UniProt ID:	P16152
Synonyms:	CBR; hCBR1; PG-9-KR; SDR21C1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)