

Product datasheet for **HP205566**

CD9 Human qPCR Primer Pair (NM_001769)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	928
Forward Sequence:	TCGCCATTGAAATAGCTGCGGC
Reverse Sequence:	CGCATAGTGGATGGCTTTCAGC
ACCN:	<u>NM_001769</u> , <u>NM_001769.1</u> , <u>NM_001769.2</u> , <u>NM_001769.3</u> , <u>BC011988</u> , <u>BC011988.1</u> , <u>BG291377</u> , <u>BM763066</u> , <u>NM_001769.4</u>
UniProt ID:	<u>P21926</u>
Synonyms:	BTCC-1; DRAP-27; MIC3; MRP-1; TSPAN-29; TSPAN29
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)