

Product datasheet for **HP205570**

CD48 Human qPCR Primer Pair (NM_001778)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	962
Forward Sequence:	GCTTGAACCACCCTTATGCCAC
Reverse Sequence:	CGTGACCACTAGCCAACCTTGCA
ACCN:	<u>BC030224</u> , <u>NM_001778</u> , <u>NM_001778.1</u> , <u>NM_001778.2</u> , <u>NM_001778.3</u> , <u>BC030224.1</u> , <u>BC016182</u> , <u>BC016182.1</u> , <u>BG505217</u> , <u>BT019813</u>
UniProt ID:	<u>P09326</u>
Synonyms:	BCM1; BLAST; BLAST1; hCD48; mCD48; MEM-102; SLAMF2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)