

Product datasheet for **HP205580**

N Cadherin (CDH2) Human qPCR Primer Pair (NM_001792)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1000
Forward Sequence:	CCTCCAGAGTTTACTGCCATGAC
Reverse Sequence:	GTAGGATCTCCGCACTGATTC
ACCN:	NM_001792 , NM_001792.1 , NM_001792.2 , NM_001792.3 , NM_001792.4 , BC036470 , BC036470.1 , BC021695 , BG197310 , BX648579 , NM_001792.5
UniProt ID:	P19022
Synonyms:	ACOGS; ARVD14; CD325; CDHN; CDw325; NCAD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)