

Product datasheet for **HP205586**

CDK7 Human qPCR Primer Pair (NM_001799)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1022
Forward Sequence:	GCACACCAACTGAGGAACAGTG
Reverse Sequence:	AAGTCGTCTCCTGCTGCACTGA
ACCN:	<u>NM_001799</u> , <u>NM_001799.1</u> , <u>NM_001799.2</u> , <u>NM_001799.3</u> , <u>BC005298</u> , <u>BC005298.1</u> , <u>BC000834</u> , <u>BI520184</u> , <u>BU660651</u> , <u>NM_001799.4</u>
UniProt ID:	<u>P50613</u>
Synonyms:	CAK; CAK1; CDKN7; HCAK; MO15; p39MO15; STK1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)