

## Product datasheet for **HP205612**

### CLPS Human qPCR Primer Pair (NM\_001832)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1208
Forward Sequence:	TTGCTGCCAGCATTCAAGTGCG
Reverse Sequence:	TCTTGTCTCCCTCACAGGTCAG
ACCN:	<a href="#">NM_001832</a> , <a href="#">NM_001832.1</a> , <a href="#">NM_001832.2</a> , <a href="#">NM_001832.3</a> , <a href="#">BC017897</a> , <a href="#">BC007061</a> , <a href="#">BC025693</a> , <a href="#">BC043489</a> , <a href="#">BI439786</a> , <a href="#">BP321569</a> , <a href="#">BP322637</a> , <a href="#">BP324927</a> , <a href="#">BQ271978</a> , <a href="#">BT006812</a> , <a href="#">NM_001832.4</a>
UniProt ID:	<a href="#">P04118</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)