

Product datasheet for **HP205623**

Collagen II (COL2A1) Human qPCR Primer Pair (NM_001844)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1280
Forward Sequence:	CCTGGCAAAGATGGTGAGACAG
Reverse Sequence:	CCTGGTTTTCCACCTTCACCTG
ACCN:	<u>BC007252</u> , <u>NM_001844</u> , <u>NM_001844.1</u> , <u>NM_001844.2</u> , <u>NM_001844.3</u> , <u>NM_001844.4</u> , <u>BC007252.1</u> , <u>BC116449</u> , <u>BT007205</u>
UniProt ID:	<u>P02458</u>
Synonyms:	ANFH; AOM; COL11A3; SEDC; STL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)