

Product datasheet for **HP205627**

Collagen VI (COL6A1) Human qPCR Primer Pair (NM_001848)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1291
Forward Sequence:	GCCTTCCTGAAGAATGTCACCG
Reverse Sequence:	TCCAGCAGGATGGTGATGTCAG
ACCN:	<u>NM_001848</u> , <u>NM_001848.1</u> , <u>NM_001848.2</u> , <u>BC052575</u> , <u>BC005159</u> , <u>BC022236</u> , <u>BC032821</u> , <u>BI769554</u> , <u>BM686857</u> , <u>BM707144</u>
UniProt ID:	<u>P12109</u>
Synonyms:	BTHLM1; OPLL; UCHMD1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)