

Product datasheet for **HP205628**

Collagen VI (COL6A2) Human qPCR Primer Pair (NM_001849)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1292
Forward Sequence:	CGTGGAGACTCAGGACAGCCA
Reverse Sequence:	CCTTTCAAGCCAAAGTCGCCTC
ACCN:	BC002484 , NM_001849 , NM_001849.1 , NM_001849.2 , NM_001849.3 , BC065509 , BC065509.1 , BM989656 , BT006913 , NM_001849.4
UniProt ID:	P12110
Synonyms:	BTHLM1; PP3610; UCMD1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)