

Product datasheet for **HP205643**

COX7A2 Human qPCR Primer Pair (NM_001865)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1347
Forward Sequence:	CGGAGAAGCAAAAAGTGTCCAG
Reverse Sequence:	CCAACTGTAAGAATCATGGTGGC
ACCN:	NM_001865 , NM_001865.1 , NM_001865.2 , NM_001865.3 , BC100852 , BC100853 , BC100854 , BC101826 , BC101828 , BC133654 , NM_001865.4
UniProt ID:	P14406
Synonyms:	COX7AL; COX7AL1; COXVIIa-L; COXVIIAL; VIIAL
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)