

## Product datasheet for **HP205651**

### Carboxypeptidase H (CPE) Human qPCR Primer Pair (NM\_001873)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1363
Forward Sequence:	TTGCTCCTGAGACCAAGGCTGT
Reverse Sequence:	GAGGAGCTGTATTCGTGAGCAC
ACCN:	<a href="#">NM_001873</a> , <a href="#">NM_001873.1</a> , <a href="#">NM_001873.2</a> , <a href="#">NM_001873.3</a> , <a href="#">BC053612</a> , <a href="#">BC053612.1</a> , <a href="#">BC033866</a> , <a href="#">BM271829</a> , <a href="#">NM_001873.4</a>
UniProt ID:	<a href="#">P16870</a>
Synonyms:	CPH; IDDHH
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)