

## Product datasheet for **HP205677**

### Cathepsin B (CTSB) Human qPCR Primer Pair (NM\_001908)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1508
Forward Sequence:	GCTTCGATGCACGGGAACAATG
Reverse Sequence:	CATTGGTGTGGATGCAGATCCG
ACCN:	<u><a href="#">NM_001908</a></u> , <u><a href="#">NM_001908.1</a></u> , <u><a href="#">NM_001908.2</a></u> , <u><a href="#">NM_001908.3</a></u> , <u><a href="#">NM_001908.4</a></u> , <u><a href="#">BC095408</a></u> , <u><a href="#">BC095408.1</a></u> , <u><a href="#">BC010240</a></u> , <u><a href="#">BG770805</a></u> , <u><a href="#">BM008741</a></u> , <u><a href="#">BX647765</a></u> , <u><a href="#">BX647952</a></u> , <u><a href="#">BX648239</a></u> , <u><a href="#">NM_001908.5</a></u>
UniProt ID:	<u><a href="#">P07858</a></u>
Synonyms:	APPS; CPSB; RECEUP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)