

## Product datasheet for **HP205681**

### Cathepsin L (CTSL) Human qPCR Primer Pair (NM\_001912)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1514
Forward Sequence:	GAAAGGCTACGTGACTCCTGTG
Reverse Sequence:	CCAGATTCTGCTCACTCAGTGAG
ACCN:	<u><a href="#">NM_001912</a></u> , <u><a href="#">NM_001912.1</a></u> , <u><a href="#">NM_001912.2</a></u> , <u><a href="#">NM_001912.3</a></u> , <u><a href="#">NM_001912.4</a></u> , <u><a href="#">BC012612</a></u> , <u><a href="#">BC012612.1</a></u> , <u><a href="#">BC142983</a></u> , <u><a href="#">BI824449</a></u> , <u><a href="#">BI829385</a></u> , <u><a href="#">BM680778</a></u> , <u><a href="#">BX537395</a></u> , <u><a href="#">BX647101</a></u> , <u><a href="#">BX647102</a></u> , <u><a href="#">BX647413</a></u> , <u><a href="#">BX647434</a></u> , <u><a href="#">BX647435</a></u> , <u><a href="#">BX648848</a></u> , <u><a href="#">BX648849</a></u> , <u><a href="#">BX649140</a></u> , <u><a href="#">NM_001912.5</a></u>
UniProt ID:	<u><a href="#">P07711</a></u>
Synonyms:	CATL; CTSL1; MEP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)