

Product datasheet for **HP205686**

DBT Human qPCR Primer Pair (NM_001918)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1629
Forward Sequence:	GCACCACTGATCTTACAGGAGG
Reverse Sequence:	CCAATGGCTACTTCAGGTGGCA
ACCN:	<u>NM_001918</u> , <u>NM_001918.1</u> , <u>NM_001918.2</u> , <u>NM_001918.3</u> , <u>BC016675</u> , <u>BC030002</u> , <u>BT007372</u> , <u>NM_001918.4</u>
UniProt ID:	<u>P11182</u>
Synonyms:	BCATE2; BCKAD-E2; BCKADE2; BCKDH-E2; BCOADC-E2; E2; E2B
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)