

Product datasheet for **HP205687**

DCI (ECI1) Human qPCR Primer Pair (NM_001919)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1632
Forward Sequence:	TGAGCCTGGAGTTTCTGACGGA
Reverse Sequence:	ACACATCTCCGTCAGGTCCAGG
ACCN:	NM_001919 , NM_001919.1 , NM_001919.2 , NM_001919.3 , BC000762 , BC000762.2 , BC002746 , BC009631 , BC019316 , BC020228 , BQ773923 , NM_001919.4
UniProt ID:	P42126
Synonyms:	DCI
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)