

Product datasheet for **HP205718**

Endothelin 2 (EDN2) Human qPCR Primer Pair (NM_001956)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1907
Forward Sequence:	CGTGTTCCAGACTGGCAAGACA
Reverse Sequence:	CCTCCTGTTGTCGCTTGGCAA
ACCN:	<u>NM_001956</u> , <u>NM_001956.1</u> , <u>NM_001956.2</u> , <u>NM_001956.3</u> , <u>NM_001956.4</u> , <u>BC034393</u> , <u>BC034393.1</u> , <u>BI761655</u> , <u>NM_001956.5</u>
UniProt ID:	<u>P20800</u>
Synonyms:	ET-2; ET2; PPET2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)