

Product datasheet for **HP205726**

EGR4 Human qPCR Primer Pair (NM_001965)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1961
Forward Sequence:	CGCCTAGGTCTGTGCGCTGC
Reverse Sequence:	CGGAAAACCTCGCTAAGGTGGAG
ACCN:	<u>NM_001965</u> , <u>NM_001965.2</u> , <u>NM_001965.3</u> , <u>BC034462</u> , <u>BC136252</u> , <u>BC143379</u> , <u>BX091371</u>
UniProt ID:	<u>Q05215</u>
Synonyms:	AT133; NGFI-C; NGFIC; PAT133
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)