

Product datasheet for **HP205730**

EIF5 Human qPCR Primer Pair (NM_001969)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1983
Forward Sequence:	AATGGCTCCGTATCCAGCAGTG
Reverse Sequence:	GCTTCCTCAGTTGTATCTTCTCC
ACCN:	<u>NM_001969</u> , <u>NM_001969.1</u> , <u>NM_001969.2</u> , <u>NM_001969.3</u> , <u>NM_001969.4</u> , <u>BC032866</u> , <u>BC000106</u> , <u>BC007728</u> , <u>BC015630</u> , <u>BC041083</u> , <u>BC107884</u> , <u>BF979642</u> , <u>BI459720</u> , <u>BX429104</u> , <u>BX537367</u> , <u>NM_001969.5</u>
UniProt ID:	<u>P55010</u>
Synonyms:	EIF-5; EIF-5A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)