

## Product datasheet for **HP205736**

### NSE (ENO2) Human qPCR Primer Pair (NM\_001975)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2026
Forward Sequence:	CTGTATCGCCACATTGCTCAGC
Reverse Sequence:	AGCTTGTTGCCAGCATGAGAGC
ACCN:	<u>NM_001975, NM_001975.1, NM_001975.2, BC002745, BC002745.2, BT007383, NM_001975.3</u>
UniProt ID:	<u>P09104</u>
Synonyms:	HEL-S-279; NSE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)