

## Product datasheet for **HP205741**

### EPS15 Human qPCR Primer Pair (NM\_001981)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Gene ID:          | 2060   |
| Forward Sequence: | AACCTGTGCCAATGTCCTTGCC   |
| Reverse Sequence: | TTCACGGACCTCCAATCCAGAC   |
| ACCN:             | <u>NM_001981, NM_001981.1, NM_001981.2, BC054006, BM681433, BX647676, NM_001981.3</u>  |
| UniProt ID:       | <u>P42566</u>  |
| Synonyms:         | AF-1P; AF1P; MLLT5   |
| Components:       | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | Store at -20°C.  |
| Stability:        | The primer mix is stable for one year from date of shipping.   |



[View online »](#)