

Product datasheet for **HP205750**

Tissue Factor (F3) Human qPCR Primer Pair (NM_001993)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2152
Forward Sequence:	CAGAGTTCACACCTTACCTGGAG
Reverse Sequence:	GTTGTTCCCTTCTGACTAAAGTCCG
ACCN:	<u>NM_001993</u> , <u>NM_001993.1</u> , <u>NM_001993.2</u> , <u>NM_001993.3</u> , <u>NM_001993.4</u> , <u>BC011029</u> , <u>BC011029.1</u> , <u>BM769013</u> , <u>BM981095</u> , <u>BT019808</u> , <u>NM_001993.5</u>
UniProt ID:	<u>P13726</u>
Synonyms:	CD142; TF; TFA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)