

Product datasheet for **HP205762**

FGF2 Human qPCR Primer Pair (NM_002006)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2247
Forward Sequence:	AGCGGCTGTACTGCAAAAACGG
Reverse Sequence:	CCTTTGATAGACACAACCTCTCTC
ACCN:	<u>NM_002006</u> , <u>NM_002006.1</u> , <u>NM_002006.2</u> , <u>NM_002006.3</u> , <u>NM_002006.4</u> , <u>BC111737</u> , <u>BC166646</u> , <u>BE439775</u> , <u>BP292299</u> , <u>BU501243</u> , <u>NM_002006.5</u>
UniProt ID:	<u>P09038</u>
Synonyms:	BFGF; FGF-2; FGFB; HBGF-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)