

Product datasheet for **HP205769**

FKBP52 (FKBP4) Human qPCR Primer Pair (NM_002014)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2288
Forward Sequence:	TGACTCCAGTCTGGATCGCAAG
Reverse Sequence:	CTGGTTTGCAGGTGATGTGGCA
ACCN:	NM_002014 , NM_002014.1 , NM_002014.2 , NM_002014.3 , BC001786 , BC001786.1 , BC002887 , BC007924 , BQ083177 , BQ446607 , BU685040 , NM_002014.4
UniProt ID:	Q02790
Synonyms:	FKBP51; FKBP52; FKBP59; HBI; Hsp56; p52; PPlase
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)