

Product datasheet for **HP205805**

GFAP Human qPCR Primer Pair (NM_002055)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2670
Forward Sequence:	CTGGAGAGGAAGATTGAGTCGC
Reverse Sequence:	ACGTCAAGCTCCACATGGACCT
ACCN:	NM_002055 , NM_002055.1 , NM_002055.2 , NM_002055.3 , NM_002055.4 , BC041765 , BC012228 , BC013596 , BC062609 , BC127871 , BF341765 , BF526347 , BM687259 , BM931499 , NM_002055.5
UniProt ID:	P14136
Synonyms:	ALXDRD
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)