

## Product datasheet for **HP205806**

### GFPT1 Human qPCR Primer Pair (NM\_002056)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2673
Forward Sequence:	CCCTCTGTTGATTGGTGTACGG
Reverse Sequence:	GGAAAAGGCAGGTTGTGCTGTC
ACCN:	<u><a href="#">NM_002056</a></u> , <u><a href="#">NM_002056.1</a></u> , <u><a href="#">NM_002056.2</a></u> , <u><a href="#">NM_002056.3</a></u> , <u><a href="#">BC045641</a></u> , <u><a href="#">BC045641.1</a></u> , <u><a href="#">BG723175</a></u> , <u><a href="#">BI560595</a></u> , <u><a href="#">BM744371</a></u> , <u><a href="#">NM_002056.4</a></u>
UniProt ID:	<u><a href="#">Q06210</a></u>
Synonyms:	CMS12; CMSTA1; GFA; GFAT; GFAT 1; GFAT1; GFAT1m; GFPT; GFPT1L; MSLG
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)