

Product datasheet for **HP205811**

GLP1R Human qPCR Primer Pair (NM_002062)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2740
Forward Sequence:	CTACGCACTCTCCTTCTCTGCT
Reverse Sequence:	CGGACAATGCTCGCAGGATGAA
ACCN:	NM_002062 , NM_002062.1 , NM_002062.2 , NM_002062.3 , NM_002062.4 , BC112126 , BC112126.1 , BC111843 , BC113493 , NM_002062.5
UniProt ID:	P43220
Synonyms:	GLP-1; GLP-1-R; GLP-1R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)