

Product datasheet for **HP205818**

GNAQ Human qPCR Primer Pair (NM_002072)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2776
Forward Sequence:	GATCAGAGCCATGGACACACTC
Reverse Sequence:	GCAGACACCTTCTCCACATCAAC
ACCN:	NM_002072 , NM_002072.1 , NM_002072.2 , NM_002072.4 , BC057777 , BC057777.1 , BC047048 , BC067850 , BC069520 , BC075096 , BC075097 , NM_002072.5
UniProt ID:	P50148
Synonyms:	CMC1; G-ALPHA-q; GAQ; SWS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)