

Product datasheet for **HP205819**

GNAZ Human qPCR Primer Pair (NM_002073)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2781
Forward Sequence:	AGCGGCTACGACCTGAAACTCT
Reverse Sequence:	GGAAGAGGATGAGTGAGGTGTTG
ACCN:	<u>NM_002073</u> , <u>NM_002073.1</u> , <u>NM_002073.2</u> , <u>NM_002073.3</u> , <u>BC037333</u> , <u>BC037333.1</u> , <u>BC026342</u> , <u>BC078163</u> , <u>BC096828</u> , <u>BU728170</u> , <u>NM_002073.4</u>
UniProt ID:	<u>P19086</u>
Synonyms:	gz-alpha
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)