

Product datasheet for **HP205824**

Golgin subfamily A member 4 (GOLGA4) Human qPCR Primer Pair (NM_002078)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2803
Forward Sequence:	AGTCGGGTGAAACAGGAGGTTG
Reverse Sequence:	GCTTCTGGTCAGTTCCTGCTC
ACCN:	<u>NM_002078</u> , <u>NM_002078.1</u> , <u>NM_002078.2</u> , <u>NM_002078.3</u> , <u>NM_002078.4</u> , <u>BC023648</u> , <u>BC045543</u> , <u>BC050444</u> , <u>BC172405</u> , <u>BX640663</u> , <u>NM_002078.5</u>
UniProt ID:	<u>Q13439</u>
Synonyms:	CRPF46; GCP2; GOLG; MU-RMS-40.18; p230
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)