

Product datasheet for **HP205836**

TFIIE beta (GTF2E2) Human qPCR Primer Pair (NM_002095)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Gene ID: | 2961 |
| Forward Sequence: | GAACATGGAGGATCGTCAGGCT |
| Reverse Sequence: | TATCTCCTCGCTGATGCCGTGT |
| ACCN: | NM_002095 , NM_002095.1 , NM_002095.2 , NM_002095.3 , NM_002095.4 , NM_002095.5 , BC030572 , BC030572.2 , BC001775 , BX421791 , NM_002095.6 |
| UniProt ID: | P29084 |
| Synonyms: | FE; TF2E2; TFIIE-B; TTD6 |
| Components: | 1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | Store at -20°C. |
| Stability: | The primer mix is stable for one year from date of shipping. |



[View online »](#)