

Product datasheet for **HP205841**

GYPB Human qPCR Primer Pair (NM_002100)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	2994
Forward Sequence:	GGCTGGTATTATTGGAACGATCC
Reverse Sequence:	GGTTCTAGGCAAGATCAGGCAG
ACCN:	NM_002100 , NM_002100.1 , NM_002100.2 , NM_002100.3 , NM_002100.4 , NM_002100.5 , BC121077 , BC069310 , BC121078 , NM_002100.6
UniProt ID:	P06028
Synonyms:	CD235b; GPB; GYP; GYPA; MNS; PAS-3; SS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)