

Product datasheet for **HP205855**

HIVEP1 Human qPCR Primer Pair (NM_002114)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3096
Forward Sequence:	ACGAACCCATACAGATGTCCGC
Reverse Sequence:	CCTACTGAGACGCCTAAATCCAC
ACCN:	<u>NM_002114, NM_002114.1, NM_002114.2, NM_002114.3, BC028111, BC140816</u>
UniProt ID:	<u>P15822</u>
Synonyms:	CIRIP; CRYBP1; GAAP; MBP-1; PRDII-BF1; Schnurri-1; ZAS1; ZNF40; ZNF40A
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)