

Product datasheet for **HP205861**

HLA DOB (HLA-DOB) Human qPCR Primer Pair (NM_002120)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3112
Forward Sequence:	CCAGATGCTGAGCAGTGGAAACA
Reverse Sequence:	GGTACACTGTCACCTCTGGTTG
ACCN:	BC020226 , NM_002120 , NM_002120.1 , NM_002120.2 , NM_002120.3 , BC006097 , BC006097.1 , BM152694 , BQ010608 , NM_002120.4
UniProt ID:	P13765
Synonyms:	DOB; HLA_DOB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)