

Product datasheet for **HP205883**

HOXB2 Human qPCR Primer Pair (NM_002145)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3212
Forward Sequence:	GAAATCCGCCAAGAAACCCAGC
Reverse Sequence:	CTGCGTGTGGTGAAGCCGTG
ACCN:	<u>NM_002145</u> , <u>NM_002145.1</u> , <u>NM_002145.2</u> , <u>NM_002145.3</u> , <u>BC074806</u> , <u>BC074806.2</u> , <u>BC038968</u> , <u>BC074805</u>
UniProt ID:	<u>P14652</u>
Synonyms:	Hox-2.8; HOX2; HOX2H; K8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)