

Product datasheet for **HP205884**

HOXB3 Human qPCR Primer Pair (NM_002146)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3213
Forward Sequence:	CCTGGATGAAAGAGTCGAGGCA
Reverse Sequence:	CCGCTGCCACCACTGCCTCC
ACCN:	<u>NM_002146, NM_002146.1, NM_002146.2, NM_002146.3, NM_002146.4, BC172242, BC172550</u>
UniProt ID:	<u>P14651</u>
Synonyms:	Hox-2.7; HOX2; HOX2G
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)