

## Product datasheet for **HP205891**

### HSD17B2 Human qPCR Primer Pair (NM\_002153)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3294
Forward Sequence:	TCCAACCTGGAGGCTTCCTAAC
Reverse Sequence:	GCTGTGCTAAGATGTAGTCTGG
ACCN:	<a href="#">NM_002153</a> , <a href="#">NM_002153.1</a> , <a href="#">NM_002153.2</a> , <a href="#">BC009581</a> , <a href="#">BC009581.1</a> , <a href="#">BC059170</a> , <a href="#">BM670647</a> , <a href="#">BP270993</a> , <a href="#">NM_002153.3</a>
UniProt ID:	<a href="#">P37059</a>
Synonyms:	EDH17B2; HSD17; SDR9C2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)