

Product datasheet for **HP205918**

IL1RAP Human qPCR Primer Pair (NM_002182)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3556
Forward Sequence:	CTGAGGATCTCAAGCGCAGCTA
Reverse Sequence:	AGCAGGACTGTGGCTCCAAAAC
ACCN:	NM_002182 , NM_002182.1 , NM_002182.2 , NM_002182.3 , BC053621 , BC053621.1 , BC016141 , BG032519 , BM194193 , BM760155 , BM973472 , NM_002182.4
UniProt ID:	Q9NPH3
Synonyms:	C3orf13; IL-1RAcP; IL1R3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)