

## Product datasheet for **HP205921**

### IL7R alpha (IL7R) Human qPCR Primer Pair (NM\_002185)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3575
Forward Sequence:	ATCGCAGCACTCACTGACCTGT
Reverse Sequence:	TCAGGCACTTTACCTCCACGAG
ACCN:	<u><a href="#">NM_002185</a></u> , <u><a href="#">NM_002185.1</a></u> , <u><a href="#">NM_002185.2</a></u> , <u><a href="#">NM_002185.3</a></u> , <u><a href="#">BC067537</a></u> , <u><a href="#">BC067537.1</a></u> , <u><a href="#">BC020717</a></u> , <u><a href="#">BC067538</a></u> , <u><a href="#">BC067539</a></u> , <u><a href="#">BC067540</a></u> , <u><a href="#">BC069999</a></u> , <u><a href="#">NM_002185.5</a></u>
UniProt ID:	<u><a href="#">P16871</a></u>
Synonyms:	CD127; CDW127; IL-7R-alpha; IL7RA; ILRA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)