

Product datasheet for **HP205933**

Aconitase 1 (ACO1) Human qPCR Primer Pair (NM_002197)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	48
Forward Sequence:	TCCTCAGGTGATTGGCTACAGG
Reverse Sequence:	TCGGTCAGCAATGGACAACCTGG
ACCN:	<u>NM_002197</u> , <u>NM_002197.1</u> , <u>NM_002197.2</u> , <u>BC018103</u> , <u>BC018103.1</u> , <u>BP287929</u> , <u>BX538349</u> , <u>NM_002197.3</u>
UniProt ID:	<u>P21399</u>
Synonyms:	ACONS; HEL60; IREB1; IREBP; IREBP1; IRP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)