

## Product datasheet for **HP205936**

### ISG20 Human qPCR Primer Pair (NM\_002201)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3669
Forward Sequence:	ACACGTCCACTGACAGGCTGTT
Reverse Sequence:	ATCTTCCACCGAGCTGTGTCCA
ACCN:	<a href="#">NM_002201</a> , <a href="#">NM_002201.1</a> , <a href="#">NM_002201.2</a> , <a href="#">NM_002201.3</a> , <a href="#">NM_002201.4</a> , <a href="#">NM_002201.5</a> , <a href="#">BC007922</a> , <a href="#">BC007922.2</a> , <a href="#">BC016341</a> , <a href="#">BG759382</a> , <a href="#">BI838264</a> , <a href="#">BM691295</a> , <a href="#">BT006952</a> , <a href="#">BX339191</a>
UniProt ID:	<a href="#">Q96AZ6</a>
Synonyms:	CD25; HEM45
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)